



The WVSM is perfect in situations where large, traditional single-patient monitors don't make sense.



WVSM is an adaptable and scalable monitoring solution for a wide variety of field applications



WVSM[®] MULTI-PATIENT MONITORING

The small and portable *Wireless Vital Signs Monitor (WVSM)* collects the data you need to make critical decisions, even in remote, dangerous and multi-patient environments.

The WVSM allows a single caregiver to monitor up to 20 patients simultaneously on a tablet, PC or smart phone.

With quick and easy attachment to a blood pressure cuff, the WVSM is ideal for patient monitoring in multiple emergency scenarios, both civilian and military

Proven by elite military medical providers in some of the world's harshest conditions...Now available to you.



FDA cleared

Different jobs require different tools. Finally there is a monitoring alternative.

- Weighs just one pound
- A fraction of the cost of a larger single-patient monitor
- Automatically records and transmits 4.5 hours of patient history
- Up to 7.5 hours of battery life
- Measures NIBP, HR, SpO2 and ECG
- Color-coded vitals provide the data medics need without overwhelming first responders

APPLICATIONS

FIRE REHAB

Patient-worn monitor allows for freedom of movement in the rehabilitation area.

EXTRICATION

Remove responders from crowded, dangerous patient vicinity.

TRIAGE/MCI

Monitors up to 20 patients on one compatible laptop, tablet or smart phone.

EMERGENCY DEPARTMENT

Monitor multiple ambulatory patients simultaneously in crowded waiting rooms.

HAZMAT

Monitor responders inside level A suits.

HIGH ANGLE RESCUE

Fits inside a stokes litter. Does not interfere with climbing harness.

INFECTIOUS DISEASE

Wireless transmission allows caregiver to monitor patient in isolation, while minimizing volumes of contaminated equipment.

SEARCH AND RESCUE

Easy to carry, ruggedized to military standards, 7.5 hours of battery life and 4.5 hours of saved data.

STRUCTURAL COLLAPSE

Small size and rugged construction goes where others can not.

Attaches to a blood pressure cuff, streamlining patient movement



WVSM FEATURES

Patient data is gathered automatically thus freeing a provider for other duties. Multiple providers can also access patient information on-scene. EMTs and Incident Commanders can access the same information in near real-time to make informed decisions about triage and patient transport. All patient data is color-coded based on user set parameters enabling easy triage.

- Highly portable
- Wired or wireless connectivity
- Red, yellow and green color coding
- Wireless range of up to 200 yards
- Turn key solution with all consumables ready to deploy
- Stays with the patient from point-of-injury until treatment

CONTACT US TODAY TO GET MORE INFORMATION OR TO SCHEDULE A DEMO.



sales@athenagtx.com



515.288.3360



athenagtx.com