

TD-1107

coretemp⁺


EAR INFRARED THERMOMETER




Apply Probe Filter



Measuring



Eject Probe Filter



Illuminated Display

 **Bluetooth[®]**
SMART



Probe filter

Hygienic and disposable probe filter for medical applications. Provided with one-hand easy operation of probe filter attachment and ejection without touching probe filter directly to reduce cross contamination.



Illuminated Display

Illuminated display with dual-colored indicator for Hi/Low and fever guidance.



Accurate and Fast

Guiding indicator to ensure proper and stable probe-tip placement to obtain clinical accuracy in seconds. The advanced measuring technology captures and analyzes multiple measurements and selects the most accurate results.



THERMOMETER SPECIFICATIONS

Measurement Range	32.0°C (89.6°F) ~43.0°C(109.4°F)
Accuracy	±0.2°C (±0.4°F) for the range of 36.0°C to 39.0°C (96.8°F to 102.2°F) ±0.3°C (±0.5°F) for the range of < 36.0°C (96.8°F) or > 39.0°C (102.2°F)
Operating Condition	10°C to 40°C (50°F to 104°F), 95% RH or less non condensing
Storage/Transportation Condition	-20°C to 60°C (-4°F to 140°F), 95% RH or less non condensing
Resolution	0.1°C / 0.1°F
Measurement Units	°C or °F Switchable
LCD Back Light	32°C 38°C 43°C Lo Normal Fever Hi
Probe Cover Detector	Equipped with Auto power-on function
Power source	AA x 2 alkaline battery
Battery life	At Least 3000 Measurements
Memory Capacity	10 Measurements
Power Saving	30 seconds
Communication	Bluetooth Smart (optional)
Dimension	
Weight	
Units per Carton	
Net Weight per Carton	
Gross Weight per Carton	
Carton Size	



PROBE COVER SPECIFICATIONS

Probe Cover Model	2300-ACC-007
Probe Cover Tray Dimension	148.76 (L) x 46 (W) x 21.7 (H) mm
Probe Cover Unit per Box	320 pcs
Probe Cover Box per Carton	6 pcs
Probe Cover Carton Size	31.5 (L) x 30 (W) x 15 (H) cm

Politek Medikal Cihazlar A.Ş.

Yenibosna Mah. M. Nevzat Pisak Cad. No:4/4
Ic Kapi No:2 Bahcelievler / Istanbul / TURKEY

Tel : +90-212-551 01 46 Fax : +90-212-551 03 11