

MAXIASPEED

OMRON



GARANTİ

3
YIL



HASTANE KULLANIMI İÇİN TEKERLEKLİ CERRAHİ ASPIRATÖR WHEELED SURGICAL ASPIRATOR FOR HOSPITAL USE

TEKNİK ÖZELLİKLER

VAKUM SEVİYESİ	Yüksek Vakum
ÇALIŞMA DÖNGÜSÜ	Sürekli
HAVA AKIŞI (Versiyon 6.xx)	60 lt/dk
HAVA AKIŞI (Versiyon 9.xx)	90 lt/dk
AYARLANABİLİR VAKUM SEVİYESİ	0 ÷ 0,90 bar (0 ÷ - 90 kPa)
BOYUTLAR	470 x 580 x 560 mm
AĞIRLIK (Versiyon 6.xx)	15 kg
AĞIRLIK (Versiyon 9.xx)	17 kg
MEDİKAL SINIFI	IIa

TECHNICAL DATA

HIGH FLOW	High vacuum
OPERATING CYCLE	Continuous
AIR FLOW (version 6xx)	60 l/min
AIR FLOW (version 9xx)	90 l/min
ADJUSTABLE VACUUM LEVEL	0 ÷ 0.90 bar (0 ÷ - 90 kPa)
DIMENSIONS	470 x 580 x 560 mm
WEIGHT (version 6xx)	15 kg
(version 9xx)	17 kg
MDD CLASS	IIa

HASTANELERDE ve AMELİYATHANELERDE VÜCUT SIVILARININ ASPIRASYONUNDA KULLANILAN, MERKEZİ ELEKTRİKLE ÇALIŞAN CERRAHİ ASPIRATÖR. TAŞIMA

KOLAYLIĞI ve SÜREKLİ KULLANIM İÇİN GELİŞTİRİLMİŞTİR.

2 ADEDİ FRENLİ 4 ANTİSTATİK TEKERLEKLE - TELESKOPIK TAŞIMA ELÇEKLERİ - EN SON AVRUPA GÜVENLİK STANDARTLARINA UYGUN YÜKSEK ISIYA DAYANIKLI ELEKTRİK İZOLASYONLU PLASTİK MALZEME - 2 ADET TAŞMA KORUMALI GÜVENLİK VALFLİ OTOKLAVLANABİLİR POLİKARBONAT KAVANOZ - ATOKSİK STERİL EDİLEBİLİR HORTUMLAR - 2 ADET %99 ANTİBAKTERİYAL HİDROFOBİK FİLTRE - BAR/PSİ VAKUMETRELİ VAKUM REGÜLATÖRÜ - KONTAMİNE OLMUŞ PARÇALARA TEMASI ENGELLEYEN BİR KAVANOZDAN DİĞERİNE GEÇİŞ SAĞLAYAN OTOMATİK AKIŞ YÖNÜ REGÜLATÖRÜ



SURGICAL ASPIRATOR WORKING WITH NETWORK ELECTRICITY TO BE USED FOR ASPIRATION OF BODY LIQUIDS IN OPERATING THEATRE AND HOSPITALS, IT HAS BEEN DEVELOPED TO BE EASY TO CARRY AND TO RUN CONTINUOUSLY

EQUIPPED WITH 4 ANTISTATIC WHEELS WITH 2 BRAKES - CARRYING TELESCOPIC HANDLE - HIGHLY HEAT RESISTANT ELECTRICALLY INSULATED PLASTIC MATERIAL IN CONFORMITY WITH THE LATEST EUROPEAN SAFETY STANDARDS - 2 AUTOCLAVABLE POLYCARBONATE JARS WITH SAFETY VALVE (OVERFLOW PROTECTION) - SET OF ATOXIC STERILIZABLE SILICONE TUBES - 99% ANTIBACTERIAL HYDROPHOBIC FILTERS (2) - VACUUM REGULATOR WITH BAR/PSI VACUOMETER - AUTOMATIC FLOW DIRECTION REGULATOR TO EASILY SWITCH FROM A JAR TO THE OTHER, AVOIDING TO MANAGE CONTAMINATED PARTS

MODEL 6.2P; 6.4P; 9.2P; 9.4P CERRAHİ ASPIRASYONUN UZAKTAN KONTROL EDİLMESİNİ SAĞLAYAN AYAK KONTROL PEDALI İLE ÇALIŞMAKTADIR, BU SAYEDE KULLANICININ ELLERİNİN SERBEST KALMASI SAĞLANIR

MODELS 6.2P, 6.4P, 9.2P, 9.4P ARE EQUIPPED WITH A FOOT-CONTROL SWITCH TO OBTAIN AN OPERATING ASPIRATION UNDER REMOTE CONTROL AND ENABLE THE OPERATOR TO KEEP HANDS FREE.



KAVANOZ 2000cc
2000cc JAR



KAVANOZ 4000cc
4000cc JAR



DEĞİŞİK UZUNLUKLARDA SİLİKON BAĞLANTI HORTUMLARI
CONNECTING SILICONE TUBE,
DIFFERENT MEASURES



AYAK KONTROL PEDALI
FOOT-CONTROL SWITCH



ANTİBAKTERİYAL FİLTRE MAXIASPEED
ANTIBACTERIAL FILTER MAXIASPEED

CE 0051
(dir. 93/42/EEC)

