



MADE IN Turkey

TECHNICAL SPECIFICATIONS

LAMP TYPE	Led / 1 pc
LIGHT INTENSITY	45.000 Lux / 30 cm 25.000 Lux / 50 cm
LIGHT SPOT DIAMETER	60 mm / 30 cm 90 mm / 50 cm
COLOR TEMPERATURE	4500° K
LAMP LIFE	50.000 hours
POWER CONSUMPTION	4 W
POWER INPUT	100-240 V, 50/60 Hz
LAMPHEAD DIAMETER	40 mm
LAMPHEAD ROTATION	360°
ARM TYPE	Flexible Steel Gooseneck (100 cm)
COLOR	Black / White
STANDART BASE	Mobile Stand w/5 Wheels
OPTIONAL MOUNTING	Table / Rail / Wall Mounted

