



MADE IN Turkey

TECHNICAL SPECIFICATIONS

LAMP TYPE	Led / 7 pcs
LIGHT INTENSITY	80.000 Lux / 30 cm 60.000 Lux / 50 cm
LIGHT INTENSITY CONTROL	Stepless Dimmer Control
LIGHT SPOT DIAMETER	90 mm / 30 cm 120 mm / 50 cm
INTENSITY ADJUSTMENT	Mechanical Dimmer
COLOR TEMPERATURE	4500° K
LAMP LIFE	50.000 hours
POWER CONSUMPTION	10 W
POWER INPUT	200-240 V, 50/60 Hz
LAMPHEAD DIAMETER	90 mm
LAMPHEAD ROTATION	360°
ARM TYPE	Flexible Steel Gooseneck (75 cm)
COLOR	Black / White / Blue / Pink
STANDART BASE	Mobile Stand w/5 Wheels
OPTIONAL MOUNTING	Table / Rail / Wall Mounted

